

Conant[®]

CONANT
2276.HK



1.60 Hi-Vex[®] MR-8 Plus

*Enhanced Impact Resistance! High Abbe Value!
The perfect balance of performance and comfort with the
aesthetics of a High Index lens.*

1.60 Hi-Vex[®]
MR-8 PLUS

HIGH IMPACT RESISTANT LENS


*Passes Ansi Z87.1-2015 High Impact Prescription Lens Qualification

1.60 Hi-Vex[®] MR-8 PLUS



High impact resistant lens



1.60 Hi-Vex[®] MR-8 Plus is an impact resistant High Index with an Abbe value of 41 providing superior optics.



COMFORT AND PERFORMANCE



* FULL Colt's Report Available upon Request : ANSI Z87.1-2015 High Impact Prescription Lens Material Qualification Conant Optical Lenses 1.60 MR-8 Plus SP HMC



SHANGHAI CONANT OPTICS CO., LTD.
 A 555 Chuan Da Rd, Shanghai, China 201299
 T (+86) 21 5859 8866/5859 5088
 F (+86) 21 5859 8686
 W <http://www.conantoptical.com.cn>

CONANT LENS INC.
 2255 Sewell Mill RD. Marietta GA. 30062
 PH-770-509-7860
 FX-770-579-9227
en.conantoptical.com.cn