

Conant[®]

CONANT
2276.HK



AC LENS

Security defense, more comfortable visual experience

www.conantoptical.com

LIGHTWEIGHT COMFORTABLE AND NOT FRAGILE

✓ 1.60 AC lens parameters

Product Name	1.60UV400	1.60UV++
Refractive Index (25°C)	n(d)=1.592; n(e)=1.595	n(d)=1.592; n(e)=1.595
Abbe Value	32	32
Specific Gravity (20°C)	1.15	1.15



CONANT 1.60 AC LENS
RESTORE A BETTER WORLD IN YOUR EYES

AC LENS

1.60 AC UV400 & UV++ S/F Single Vision Lens

Diam.	Base	Center Thickness	Edge Thickness	Diam.	Base	Center Thickness	Edge Thickness
70mm	PLANO	7.00	15.00	75mm	PLANO	6.00	15.00
	50	7.50	15.00		50	6.50	15.00
	100	9.00	15.00		100	7.50	15.00
	200	9.50	15.00		200	9.00	15.00
	300	8.50	13.00		300	8.00	13.00
	400	9.50	13.00		400	9.20	13.00
	500	10.50	13.00		500	10.50	13.00
	600	7.00	7.00		600	8.00	7.50
	700	8.20	7.00		700	9.10	7.50
	800	10.00	7.00		800	11.10	7.50
900	11.50	7.00	900	12.70	7.50		
1000	13.00	7.00	1000	15.10	7.50		



[AC LENS]

AC Lens Can Reduce The Blue Light Emitted By Displays And The Ultraviolet Rays Generated By Sunlight!



• COMPUTERS



• MOBILE PHONES



• TELEVISIONS



• SUN



HIGH TRANSPARENCY
TAKE CARE OF YOUR EYES
MORE AT EASE